

# Technical Datasheet

---

**Article no.: 2657**

**Acetonitril gradient for HPLC (min. 99.9 %)**

CH<sub>3</sub>CN

For laboratory use.

Parameter	Value
Water (KF)	max. 150 ppm
Non-volatile substances	max. 0.0005 %
Acidity or Alkalinity	max. 0.0005 meq/g
UV Transmission at 195 nm	min. 70 %
UV Transmission at 200 nm	min. 80 %
UV Transmission at 210 nm	min. 88 %
UV Transmission at 220 nm	min. 94 %
UV Transmission at 230 nm	min. 98 %
UV Transmission at 240 nm	min. 99 %
Colour (Hazen)	max. 10
filtered through 0.2 µm	

Version No. 02/2024-09-16