

Ordering Information

Acrodisc® Syringe Filters With PVDF Membrane, 13 mm

Part Number	Description	Pkg
4450	0.2 µm, minispikes outlet	100/pkg, 300/cs
4544	0.2 µm, minispikes outlet	1000/pkg
4452	0.45 µm, minispikes outlet	100/pkg, 300/cs
4545	0.45 µm, minispikes outlet	1000/pkg
4455	0.2 µm, male slip luer outlet	100/pkg, 300/cs
4457	0.45 µm, male slip luer outlet	100/pkg, 300/cs

Acrodisc Syringe Filters With PVDF Membrane, 25 mm

Part Number	Description	Pkg
4406	0.2 µm	50/pkg, 200/cs
4520	0.2 µm	1000/pkg
4519	0.45 µm, AutoPack™ tubes	25/pkg, 200/cs
4408	0.45 µm	50/pkg, 200/cs
4500	0.45 µm	1000/pkg

Related Products

47 mm Filter Funnels, Glass	206
AcroPrep™ 24 Filtration System	201
HPLC Mobile Phase Filtration Membranes	177
SolVac® Filter Holder	182
Stainless Steel Forceps	230, 280
Vacuum/Pressure Pumps	279



Zymark® and SOTAX® Automation Certified Acrodisc PSF Syringe Filters With PVDF Membrane, 25 mm

Part Number	Description	Pkg
AP-4795	0.2 µm, AutoPack tubes	25/pkg, 200/cs
AP-4796	0.2 µm	50/pkg, 200/cs
AP-4797	0.2 µm	1000/pkg
AP-4519	0.45 µm, AutoPack tubes	25/pkg, 200/cs
AP-4408	0.45 µm	50/pkg, 200/cs
AP-4500	0.45 µm	1000/pkg



Zymark and SOTAX Automation Certified Acrodisc PSF GxF Syringe Filters With PVDF Membrane, 25 mm

Part Number	Description	Pkg
AP-4792	GxF/0.2 µm, AutoPack tubes	25/pkg, 200/cs
AP-4793	GxF/0.2 µm	50/pkg, 200/cs
AP-4794	GxF/0.2 µm	1000/pkg
AP-4309	GxF/0.45 µm, AutoPack tubes	25/pkg, 200/cs
AP-4310	GxF/0.45 µm	50/pkg, 200/cs
AP-4308	GxF/0.45 µm	1000/pkg