

## FICHE TECHNIQUE

<b>Référence</b>	28-4065-53
<b>Désignation</b>	DNTP SET - A C G T 4x500µmoles



### Description :

- ✓ Parameter Amersham dNTP Set (100 mM each A,C,G,T)
- ✓ Order Information Custom formulation, testing, packaging, and labeling are available.
- ✓ Purity > 99% triphosphate
- ✓ Physical State Liquid
- ✓ Concentration 100 mM
- ✓ Color Colorless
- ✓ Kit contents 2'-Deoxyadenosine 5'-Triphosphate Solution; 2'-Deoxycytidine 5'-Triphosphate Solution; 2'-Deoxyguanosine 5'-Triphosphate Solution; 2'-Deoxythymidine 5'-Triphosphate Solution
- ✓ Licensing For licensing information, see the Licensing Statements page in the About Us section.
- ✓ Pack size 4 x 500 µmol

### Unité de vente :

4x500 µmol