



Amersham Solution dNTPs

Amersham DNA Polymerization Mix dNTP (A,C,G,T)

28406558

Our dNTPs offer the highest purity deoxynucleotides available for DNA amplification, dideoxy sequencing, labelling, mutagenesis, complementary DNA [\(cDNA\) synthesis](#), and protein expression profiling.

- Free from DNase, RNase, and nicking enzyme activity.
- Buffer free and ready-to-use solutions in multiple formats.
- Custom formulation, testing, packaging, and labeling are available.

Functionally tested for long Polymerase Chain Reaction (PCR) and sequencing, with greater than 99% triphosphate purity (by HPLC) for optimum performance and consistency. A [Certificate](#) of Analysis for each lot is available on request

Concentration	20 mM
Order Information	Custom formulation, testing, packaging, and labeling are available.
Purity	> 99% triphosphate
Physical State	Liquid
Color	Colorless
Licensing	For licensing information, see the Licensing Statements page in the About Us section.
Pack size	4 x 10 μ mol