

---

## His SpinTrap Kit

His SpinTrap is a prepacked, single-use spin column for purifying histidine-tagged proteins by immobilized metal ion affinity chromatography (IMAC)

His SpinTrap is a prepacked, single-use spin column for purifying histidine-tagged proteins by immobilized metal ion affinity chromatography (IMAC)

- For simple, small-scale purification of histidine-tagged\* proteins and rapid expression screening.
- Ni Sepharose High Performance in prepacked spin columns.
- Yields up to 750 µg pure protein per column.
- Allows direct purification of both clarified and unclarified cell lysates.

### His SpinTrap Kit

\* See the section for licensing information in our terms and conditions.

His SpinTrap Kit



<b>Parameter</b>	<b>His SpinTrap Kit</b>
Chromatography technique	Histidine-tagged protein purification
Binding Capacity/Column <sup>1</sup>	At least 40 mg histidine-tagged protein/ml medium
Includes	1 pack His SpinTrap and 1 pack His Buffer Kit.
Licensing	For licensing information, see the Licensing Statements page in the About Us section.
Order Information	Includes 1 pack His SpinTrap and 1 pack His Buffer Kit
Storage	4 to 30°C, 0.15% Kathon CG
Pack size	1

<sup>1</sup>Protein binding capacity is protein-to-protein dependent.