



PRODUCT CODE: 283098

Silver Sulfate solution 10 g/l in sulfuric acid for volumetric analysis
for determination of COD according to Norme Française NFT 90-101

$\text{Ag}_2\text{O}_4\text{S}$

Ag_2SO_4

M.= 311,83

TARIC 3822 00 00 00

PHYSICAL DATA: Miscible with water Thermal decom. > 340 °C • D 20/4 1,84 •

BIBLIOGRAPHY:

HAZARDOUS: RTECS: VW 4725000 •



H: H314 • H290 • H410 •

P: P305+P351+P338 • P310 • P280 •

TRANSPORT REGULATIONS: UN: 1830 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,84 kg 1kg-0,54 l

SPECIFICATIONS:

Assay (as Ag_2SO_4) 9,8 - 10,2 g/l

Suitability: for COD determination passes test