

Brand: Honeywell Fluka
Product: 31642

Hydrogen peroxide 30 % by weight H₂O₂, puriss. p.a., Reag. ISO, Reag. Ph. Eur., stabilized

Parameter	Specification	Units
assay	min. 30	%
non-volatile matter	max. 0.005	%
residue on ignition	max. 0.002	%
free acid (as H ₂ SO ₄)	max. 0.005	%
acidity	complying	
aluminium (Al)	max. 0.00005	%
arsenic (As)	max. 0.000001	%
barium (Ba)	max. 0.000005	%
beryllium (Be)	max. 0.000001	%
bismuth (Bi)	max. 0.00001	%
calcium (Ca)	max. 0.00002	%
cadmium (Cd)	max. 0.000001	%
cobalt (Co)	max. 0.000001	%
chromium (Cr)	max. 0.000002	%
copper (Cu)	max. 0.000001	%
iron (Fe)	max. 0.000005	%
germanium (Ge)	max. 0.000005	%
potassium (K)	max. 0.00001	%
lithium (Li)	max. 0.000001	%
magnesium (Mg)	max. 0.000005	%
manganese (Mn)	max. 0.000001	%
molybdenum (Mo)	max. 0.000002	%
sodium (Na)	max. 0.00001	%
nickel (Ni)	max. 0.000002	%
lead (Pb)	max. 0.000001	%
strontium (Sr)	max. 0.000001	%

Hydrogen peroxide 30 % by weight H₂O₂, puriss. p.a., Reag. ISO, Reag. Ph. Eur., stabilized

Parameter	Specification	Units
titanium (Ti)	max. 0.00001	%
thallium (Tl)	max. 0.000005	%
vanadium (V)	max. 0.000001	%
zinc (Zn)	max. 0.000005	%
zirconium (Zr)	max. 0.00001	%
chloride (Cl)	max. 0.00005	%
phosphate (PO ₄)	max. 0.0001	%
sulfate (SO ₄)	max. 0.0001	%
total N	max. 0.0002	%
organic stabilisators	max. 0.05	%