

Stomacher® Bags - Strainer Bags

Product Codes by Sampling Volume

5-80ml:: **BA6040/STR, BA6040/STR/CLR, BA6040/STR/DBL**
 80-400ml: **BA6041/STR, BA6041/LSTR, BA6141/STR**


TECHNICAL SPECIFICATIONS

- Stomacher® blender bag, with addition of strainer element
- Outer bag – Virgin grade, Polyolefin complex. Flexible and transparent
- Strainer bags – Virgin grade, microfiltered Polyolefin complex. Flexible and opaque
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenisation of the sample
- Compatible with any paddle lab blender
- Suitable for food contact. Not designed to be used for food packaging
- Irradiated sterile in sealed sachets, with batch traceable certificate
- Storage: in a temperate and dry place, away from humidity and light
- Available in 80ml and 400ml

Stomacher® 80 Biomaster Strainer Bag


P/N. **BA6040/STR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)




P/N. **BA6040/STR/CLR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)
- Additional closure tab, consisting of wire reinforced tape



P/N. **BA6040/STR/DBL**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)



Seward Ltd

Technology Centre
 Easting Close
 Worthing, BN14 8HQ
 United Kingdom
 info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
 Bohemia
 NY 11716
 USA
 info@sewardusa.com


FermionX GmbH

Münstertäler Str. 12,
 Eschbach
 DE-79427
 Germany
 info@seward.co.uk

Stomacher® 400 Strainer Bags


P/N. BA6041/STR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Straight bottom seal
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)




P/N. BA6041/LSTR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Lateral side strainer arrangement
- Interwoven microfiber filter
- Straight bottom seal
- Pack size: 500 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



P/N. BA6141/STR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Round bottom seal
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



Stomacher® is registered trademark of Seward Ltd.

SEWARD is a part of the FermionX Group. FermionX Ltd reserves the right to change and improve specifications of its products without notice. Please contact info@seward.co.uk or visit www.seward.co.uk for more information.

Seward Ltd

Technology Centre
Easting Close
Worthing, BN14 8HQ
United Kingdom
info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
Bohemia
NY 11716
USA
info@sewardusa.com

FermionX GmbH

Münstertäler Str. 12,
Eschbach
DE-79427
Germany
info@seward.co.uk