

TECHNICAL DATA SHEET

Part Number	32654-25G-HON
Description	Methyl Red

Technical specifications :

- ✓ methyl red, reag. ph. eur., acid-base indicator 25g
- ✓ CAS N°493-52-7
- ✓ Molecular weight 269.30 g/mol
- ✓ Molecular formula $(\text{CH}_3)_2\text{NC}_6\text{H}_4\text{N}=\text{NC}_6\text{H}_4\text{CO}_2\text{H}$
- ✓ Density 0.989 g/cm³
- ✓ Flash point no data available
- ✓ Material form solid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 25g

- ✓ P Statements P280

Sales unit :

25g Methyl Red
