

TECHNICAL DATA SHEET

Part Number	32661-25G-HON
Description	Phenol Red

Technical specifications :

- ✓ phenol red, indicator, acs reagent 25g
- ✓ CAS N° 143-74-8
- ✓ Molecular weight 354.38 g/mol
- ✓ Molecular formula C₁₉H₁₄O₅S
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form powder

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 25g

- ✓ P Statements P280

Sales unit :

25g Phenol Red
