

TECHNICAL DATA SHEET

Part Number	34607-30ML-HON
Description	Phenolphthalein solution

Technical specifications :

- ✓ phenolphthalein solution, indicator, reag. ph. eur., 1% in ethanol 30ml
- ✓ Molecular formula C₂₀H₁₄O₄
- ✓ no data available
- ✓ Flash point 14 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 30mL

- ✓ GHS pictograms GHS07 ; GHS08 ; GHS02
- ✓ H Statements H350 ; H341 ; H319 ; H225
- ✓ P Statements P280 ; P201 ; P210 ; P243

Sales unit :

30mL Phenolphthalein solution
