


Specification

Proteinase K, lyophilized

A3830

Solubility:	soluble (H ₂ O)
Product Code:	A3830
Product Name:	Proteinase K, lyophilized
Short Description:	delivery form: lyophilized
Specifications:	DNases/RNases: not detectable Activity: min. 30 mAnsonU/mg Appearance: white pH (1 %; H ₂ O; 20°C): 6.2 - 6.8
Hazard pictograms	
WGK:	1
Storage:	2 - 8°C
Signal Word:	Danger
GHS Symbols:	GHS07 GHS08
H Phrases:	H315 H319 H334 H335
P Phrases:	P284 P305+P351+P338 P342+P311 P403+P233 P405


AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340

Specification

Proteinase K, lyophilized

A3830

	P501
Origin:	from <i>Tritirachium album</i>
M:	27000 g/mol
CAS:	39450-01-6
EINECS:	254-457-8
CS:	35079090
Index Nr.:	647-014-00-9
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340