

# Honeywell

## CERTIFICATE OF ANALYSIS

Honeywell Specialty Chemicals Seelze GmbH  
Wunstorferstrasse 40  
Seelze, 30926  
Telefon: +49 5137 999-0  
[lab.honeywell.com](http://lab.honeywell.com)

Brand: Honeywell Riedel-de Haën  
Product: 34858  
Lot No: P110M  
Production date: 20.Apr.2025  
Expiry Date: 07.Oct.2027

Ethyl acetate CHROMASOLV, for HPLC, 99.7%

Parameter	Specification	Units	Result	Units
assay (GC)	min. 99.7	%	99.98	%
non-volatile matter	max. 0.0005	%	0.00005	%
water (Karl Fischer)	max. 0.03	%	0.002	%
free acid (as CH <sub>3</sub> COOH)	max. 0.005	%	<0.005	%
absorbance at 400 nm	max. 0.01		0.001	
absorbance at 300 nm	max. 0.01		0.001	
absorbance at 270 nm	max. 0.10		0.003	
absorbance at 260 nm	max. 0.30		0.088	
absorbance at 255 nm	max. 0.70		0.500	



filtered via 0,2µm pore size filter

QC Release Date: 02.May.2025

The minimum shelf life is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.

We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.

Particular properties of the products or the suitability for a particular area of application are not assured.

We guarantee a proper quality within our General Conditions of Sales.

Honeywell

Quality Management

This document was produced electronically and is valid without a signature