



Specifications sheet

105000-01

5L PET PharmaTainer

Doc #

LUXSPE-37

Rev #

01

Approval date

05-juin-25

5L PET PharmaTainer

Catalogue number - 105000-01

Bottle part #	105000	USP Class VI	ADCF	DMF #
Bottle material	SaniSure resin code 001	Yes	Yes	15205
Nominal volume	5.000ml			
Brim Capacity	6.000ml		Tolerance	
Bottle dimensions w/o caps	H322 x W168 x 168 mm		+/- 1.8mm	
Bottle weight w/o caps	362.2		+/- 3%	
Bottle thread O.D	63.77mm		+/- 0.3mm	
Bottle neck I.D	54.5mm		+/- 0.3mm	

Graduation	500ml - 5.000ml	500ml intervals	Tolerance	
			+/- 5%	
Sterility	C60 y irradiated X-Ray	SAL10-6	Internal Surface	
			193557mm ²	
Endotoxin	No detectable Endotoxin	Certified by lot	Internal Volume	
			5718391mm ³	
Bottle decoration	Graduation scale	Lot number	Expiration date	Data Matrix
Bottle colour	Clear			

Traceability	Lot Number	Data Matrix	Serial Number
---------------------	------------	-------------	---------------

Closure part #	300070C1-99	Closure colour	Natural	
Closure Material	SaniSure resin code 003	USP Class VI	ADCF	DMF #
		Yes	Yes	1646
Closure weight	23.6g	Tolerance		
		+/- 2g		
Closure Dimensions	H 28.3mm O.D 70.0mm	Tolerance		
		+/- 3mm		

Packaging details	Carton	Level 3 - Carton liner	Level 2	Level 1 - Sterile Barrier
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µm LDPE monolayer film Heat sealed bag	120µm PE/PA/PE multilayer film vacuum bag heat sealed

Carton dimensions	L 590 – W 400 – H 350 mm	Carton weight	Units per carton	Units per sterile barrier
		4.00kg	6	1 (Single Pack)

Palettisation	Dimensions L x W x H	Weight (inc. Pallet)	Carton per pallet	Density
Industry pallet	1200x1000x1910mm	116kg	25	0.04