

 <b>Specifications sheet</b>	<h1>202000-05</h1> <h2>2L PC PharmaTainer</h2>	<b>Doc #</b>	<b>Rev #</b>
		LUXSPE-36	3
		<b>Approval date</b>	
		30 Oct 2025	

<b>2L PC PharmaTainer</b>	<b>Catalogue number - 202000-05</b>
---------------------------	-------------------------------------

<b>Bottle part #</b>	202000	<b>USP Class VI</b>	<b>ADCF</b>	<b>DMF #</b>		
<b>Bottle material</b>	SaniSure resin code 002	Yes	Yes	15205		
<b>Nominal volume</b>	2L					
<b>Brim Capacity</b>	2.270L					
<b>Bottle dimensions w/o caps</b>	H240 - W120 - D120				<b>Tolerance</b>	+/- 1,5mm
<b>Bottle weight w/o caps</b>	202g					+/- 3 %
<b>Bottle thread O.D</b>	44mm					+/- 0,3mm
<b>Bottle neck I.D</b>	37mm		+/- 0,3mm			

<b>Graduation</b>	200ml - 2000ml	100ml Intervals	<b>Tolerance</b>	
			+/- 5%	
<b>Sterility</b>	C60 y irradiated XRAY	SAL 10-6	<b>Internal Surface</b>	
			103035 mm <sup>2</sup>	
<b>Endotoxin</b>	No detectable endotoxin	Certified by Lot	<b>Internal Volume</b>	
			2297685 mm <sup>3</sup>	
<b>Bottle decoration</b>	Graduation scale	Lot number	Expiration date	Data Matrix
<b>Bottle colour</b>	Clear			

<b>Traceability</b>	Lot Number	Data Matrix	Serial Number
---------------------	------------	-------------	---------------

<b>Closure part #</b>	300048C-99	<b>Closure colour</b>	Natural	
<b>Closure Material</b>	SaniSure resin code 003	<b>USP Class VI</b>	<b>ADCF</b>	<b>DMF #</b>
		Yes	Yes	1646
<b>Closure weight</b>	12,6g	<b>Tolerance</b>		
		+/- 2g		
<b>Closure Dimensions</b>	H 25mm - O.D 49mm	<b>Tolerance</b>		
		+/- 5%		

<b>Packaging details</b>	<b>Carton</b>	<b>Level 3 - Carton liner</b>	<b>Level 2</b>	<b>Level 1 - Sterile Barrier</b>
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µm LDPE monolayer film Heat sealed bag	120µm PE/PA multilayer film vacuum bag heat sealed

<b>Carton dimensions</b>	L 630 - W 515 - H 300 mm	<b>Carton weight</b>	<b>Units per carton</b>	<b>Units per sterile barrier</b>
		6.2kg	20	5

<b>Palettisation</b>	<b>Dimensions L x W x H</b>	<b>Weight (inc. Pallet)</b>	<b>Carton per pallet</b>	<b>Density</b>
<b>Industry pallet</b>	1200x1000x1920mm	112kg	16	0.08