

FICHE TECHNIQUE

Référence	354597
Désignation	Corning® BioCoat® FluoroBlok™



Description :

- ✓ Corning® BioCoat® FluoroBlok™ Fibronectin Inserts with 3.0 µm PET Membrane in two 24-well Plates
- ✓ Brand Corning®
- ✓ Cell Growth Area 4,2 cm²
- ✓ Membrane Material Polyethylene terephthalate (PET)
- ✓ Membrane Pore Size 3 µm
- ✓ Membrane Diameter 6,4 mm
- ✓ Packaging Format 12 Inserts/Plate/Pk, 2 Plates/Cs

Unité de vente :

24 unités
