




# Polycarbonate Culture Tubes

## DATA SHEET



Catalog #	TPC20125	TPC16125	TPC15103	TPC13100
<b>Size</b>	20x125 mm	16x125 mm	15x103 mm	13x100 mm
<b>Total Volume</b>	26 cc	15 cc	11 cc	8 cc
<b>A • Diameter</b>	20 mm	16 mm	15 mm	13 mm
<b>B • Height (w/o Cap)</b>	125 mm	125 mm	103.5 mm	100 mm
<b>C • Height (w/ Cap)</b>	127.1 mm	127.1 mm	105.6 mm	102.1 mm



TriForest Lexan Polycarbonate (PC) culture tubes are optically clear, durable and a safer alternative to glass. They are ideal for production of terminally sterilized products. The tubes are flat bottom and ideal for large scale production of culture media using automated or manual filling systems. PC tubes offer a wide range of chemical resistance as well as excellent strength for repeated use. Leak-proof caps provide security during handling and transportation. TriForest culture tubes maintain excellent pH control. They are provided non sterile with caps separately packed.

TriForest Enterprises, Inc.

6 Chrysler, Irvine, CA 92618, USA | (888) 777-9956 | (949) 380-8866 | Fax: (949) 380-9955 | [www.triforest.com](http://www.triforest.com) | [info@triforest.com](mailto:info@triforest.com)

© Copyright of TriForest Enterprises, Inc. Vol. 1

# Polycarbonate Culture Tubes

## TUBE & CAP SPECIFICATIONS

	Catalog #	TPC20125	TPC16125	TPC15103	TPC13100
Tube Specification	Transparent	YES			
	Leakproof	YES			
	Tubes Shrink Wrapped	YES			
	Particle Free	YES			
	Sterile (Gamma Irradiation)	NO			
	Temperature Range	-125° C to 125° C			

Sterilization	Autoclavable	YES			
	Maximum Autoclave Temperature	125° C			
	Autoclave Cycles	≤ 10			

	Size	20-415	13-415	13-415	13-415
Cap Specification	Standard Color	Blue			
	Cap Type	Standard Seal Cap			
	Material	Polypropylene			

Packaging	Number of Tubes / Tray	150	200	286	250
	Number of Trays / Case	6	6	6	6
	Total Tubes / Case	900	1200	1716	1500
	Case Dimension (L x W x H inches)	15.5 x 13.5 x 19.5	13.5 x 12.5 x 18.5	16.5 x 12.5 x 17	14 x 10.5 x 17.5
	Case Dimension (L x W x H cm)	39 x 34 x 50	34 x 32 x 47	42 x 32 x 43	36 x 27 x 44
	Case Gross Weight	31 lb / 14 kg	40 lb / 18 kg	34 lb / 15.5 kg	24 lb / 11 kg