



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

Nome e indirizzo
Name and address of the firm

BIOSIGMA SpA
Via Valletta, 6 – Zona Ind. Cantarana
30010 - Cona (VE) - ITALY
Tel. +39.0426.302224 - Fax +39.0426.302228
E-mail: info@biosigmaeu.com - <http://www.biosigma.com>

Dichiara sotto la propria responsabilità che i dispositivi medico-diagnostici in vitro:
Declare under its own responsibility that in vitro diagnostic medical devices:

PRODOTTI PER CRIOGENIA STERILI STERILE CRYOGENIC PRODUCTS Vedi elenco codici allegato See the enclosed codes list	Codice CND	W050301020102	Provette senza additivi in materiale plastico per analisi
	EDMA code:	51091001	Other containers for sample of human origin
	GMDN code:	47775	General specimen container IVD, no additive/medium

Classificazione secondo la direttiva relativa ai dispositivi medico-diagnostici in vitro 98/79/CE e ss.mm.ii.
Classification according to the directive on in vitro diagnostic medical devices 98/79/EC and subsequent modifications and integrations.

- Dispositivo dell'elenco A, allegato II / Device of List A, Annex II
- Dispositivo dell'elenco B, allegato II / Device of List B, Annex II
- Dispositivo per test autodiagnostico non elencato nell'allegato II / Device for self-testing not listed in Annex II
- Altro dispositivo / Other device**

Sono conformi ai requisiti essenziali di cui all'allegato I della direttiva 98/79/CE del 27 ottobre 1998 recepita in Italia con il D.Lgs 332 del 08/09/2000 e ss.mm.ii.
Are manufactured in compliance with essential requirements of Annex I of the 98/79/CE Directive dated 27th October 1998 and subsequent modifications and integrations.

I Dispositivi di cui alla lista allegata non rientrano nell'elenco A o B di cui all'Allegato II della Direttiva 98/79/CE e ss.mm.ii.
The devices as per enclosed list do not fall under list A or B of annex II of the Directive 98/79/EC and subsequent modifications and integrations.

Procedimento di valutazione della conformità: **Allegato III (escluso punto 6)**
Conformity assessment procedure: Annex III (excluding point 6)

Organo incaricato della valutazione della conformità (se consultato): **N.A.**
Notified Body (if consulted):

CONA 27/01/2026
Luogo, data / Place, date

Nome e funzione / Name and function
MARCOLIN MARTINO
Consigliere Delegato/Managing Director



La presente Dichiarazione di Conformità CE resta valida fino al 31/12/2029 in conformità al REGOLAMENTO (UE) 2024/1860 DEL PARLAMENTO EUROPEO E DEL CONSIGLIO che modifica il regolamento (UE) 2017/746 per quanto riguarda le disposizioni transitorie per i dispositivi della classe A immessi sul mercato in condizioni di sterilità.
This EC Declaration of Conformity remains valid until 31/12/2029 in accordance with REGULATION (EU) 2024/1860 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL amending regulation (EU) 2017/746 as regards the transitional provisions for devices of class A placed on the market under sterile conditions.



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
| E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
| R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

LISTA CODICI / CODES LIST

REF	DESCRIZIONE	DESCRIPTION
CL1ARBEPSTS	Provetta CryoGen® da 1,2 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BAZ	Provetta CryoGen® da 1,2 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BBI	Provetta CryoGen® da 1,2 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BG	Provetta CryoGen® da 1,2 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BRO	Provetta CryoGen® da 1,2 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BV	Provetta CryoGen® da 1,2 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BVI	Provetta CryoGen® da 1,2 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/AZ	Provetta CryoGen® da 1,2 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/BI	Provetta CryoGen® da 1,2 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/G	Provetta CryoGen® da 1,2 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/RO	Provetta CryoGen® da 1,2 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/V	Provetta CryoGen® da 1,2 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/VI	Provetta CryoGen® da 1,2 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS	Provetta CryoGen® da 1,2 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BAZ	Provetta CryoGen® da 1,2 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BBI	Provetta CryoGen® da 1,2 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>



Biosigma S.p.A.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL1ARBIPSTS/BG	Provetta CryoGen® da 1,2 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BRO	Provetta CryoGen® da 1,2 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BV	Provetta CryoGen® da 1,2 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BVI	Provetta CryoGen® da 1,2 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/AZ	Provetta CryoGen® da 1,2 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/BI	Provetta CryoGen® da 1,2 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/G	Provetta CryoGen® da 1,2 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/RO	Provetta CryoGen® da 1,2 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/V	Provetta CryoGen® da 1,2 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/VI	Provetta CryoGen® da 1,2 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEP2D	Provetta CryoGen® da 2 ml con base, tappo esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL2ARBEP2D/B	Provetta CryoGen® da 2 ml con base, tappo esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL2ARBEPSTS	Provetta CryoGen® da 2 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/B	Provetta CryoGen® da 2 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BAZ	Provetta CryoGen® da 2 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BBI	Provetta CryoGen® da 2 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BG	Provetta CryoGen® da 2 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BRO	Provetta CryoGen® da 2 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BV	Provetta CryoGen® da 2 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>



Biosigma S.p.A. a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL2ARBEPSTS/BVI	Provetta CryoGen® da 2 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/AZ	Provetta CryoGen® da 2 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/BI	Provetta CryoGen® da 2 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/G	Provetta CryoGen® da 2 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/RO	Provetta CryoGen® da 2 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/V	Provetta CryoGen® da 2 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/VI	Provetta CryoGen® da 2 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPS2D	Provetta CryoGen® da 2 ml con base, tappo interno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL2ARBIPS2D/B	Provetta CryoGen® da 2 ml con base, tappo interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS	Provetta CryoGen® da 2 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/B	Provetta CryoGen® da 2 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BAZ	Provetta CryoGen® da 2 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BI	Provetta CryoGen® da 2 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BG	Provetta CryoGen® da 2 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BRO	Provetta CryoGen® da 2 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BV	Provetta CryoGen® da 2 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BVI	Provetta CryoGen® da 2 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/AZ	Provetta CryoGen® da 2 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/BI	Provetta CryoGen® da 2 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/G	Provetta CryoGen® da 2 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL2ARBIPSTS/RO	Provetta CryoGen® da 2 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/V	Provetta CryoGen® da 2 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/VI	Provetta CryoGen® da 2 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS	Provetta CryoGen® da 2 ml, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/B	Provetta CryoGen® da 2 ml, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BAZ	Provetta CryoGen® da 2 ml, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BBI	Provetta CryoGen® da 2 ml, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BG	Provetta CryoGen® da 2 ml, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BRO	Provetta CryoGen® da 2 ml, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BV	Provetta CryoGen® da 2 ml, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BVI	Provetta CryoGen® da 2 ml, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/AZ	Provetta CryoGen® da 2 ml, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/BI	Provetta CryoGen® da 2 ml, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/G	Provetta CryoGen® da 2 ml, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/RO	Provetta CryoGen® da 2 ml, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/V	Provetta CryoGen® da 2 ml, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/VI	Provetta CryoGen® da 2 ml, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS	Provetta CryoGen® da 2 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/B	Provetta CryoGen® da 2 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BAZ	Provetta CryoGen® da 2 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BBI	Provetta CryoGen® da 2 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BG	Provetta CryoGen® da 2 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BRO	Provetta CryoGen® da 2 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BV	Provetta CryoGen® da 2 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>

Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZ02A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626





DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL2ARIPSTS/BVI	Provetta CryoGen® da 2 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®- IN BULK</i>
CL2ARIPSTS/AZ	Provetta CryoGen® da 2 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/BI	Provetta CryoGen® da 2 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/G	Provetta CryoGen® da 2 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/RO	Provetta CryoGen® da 2 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/V	Provetta CryoGen® da 2 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/VI	Provetta CryoGen® da 2 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSSN	Provetta CryoGen® da 3 ml con base, tappo a vite esterno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® - con bar-code a sequenza personalizzata	<i>CryoGen® Tube 3ml self standing, external bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>
CL3ARBEPSTS	Provetta CryoGen® da 3 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/B	Provetta CryoGen® da 3 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BAZ	Provetta CryoGen® da 3 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BBI	Provetta CryoGen® da 3 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BG	Provetta CryoGen® da 3 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BRO	Provetta CryoGen® da 3 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BV	Provetta CryoGen® da 3 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BVI	Provetta CryoGen® da 3 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/AZ	Provetta CryoGen® da 3 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/BI	Provetta CryoGen® da 3 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/G	Provetta CryoGen® da 3 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/RO	Provetta CryoGen® da 3 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/V	Provetta CryoGen® da 3 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>

Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.A. IT11080P00002626





DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL3ARBEPSTS/VI	Provetta CryoGen® da 3 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS	Provetta CryoGen® da 4 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/B	Provetta CryoGen® da 4 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BAZ	Provetta CryoGen® da 4 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BBI	Provetta CryoGen® da 4 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BG	Provetta CryoGen® da 4 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BRO	Provetta CryoGen® da 4 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BV	Provetta CryoGen® da 4 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BVI	Provetta CryoGen® da 4 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/AZ	Provetta CryoGen® da 4 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/BI	Provetta CryoGen® da 4 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/G	Provetta CryoGen® da 4 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/RO	Provetta CryoGen® da 4 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/V	Provetta CryoGen® da 4 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/VI	Provetta CryoGen® da 4 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS	Provetta CryoGen® da 4 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/B	Provetta CryoGen® da 4 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BAZ	Provetta CryoGen® da 4 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BBI	Provetta CryoGen® da 4 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BG	Provetta CryoGen® da 4 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZ02A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL4ARBIPSTS/BRO	Provetta CryoGen® da 4 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BV	Provetta CryoGen® da 4 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BVI	Provetta CryoGen® da 4 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/AZ	Provetta CryoGen® da 4 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/BI	Provetta CryoGen® da 4 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/G	Provetta CryoGen® da 4 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/RO	Provetta CryoGen® da 4 ml con base, tappo interno RED, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/V	Provetta CryoGen® da 4 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/VI	Provetta CryoGen® da 4 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS	Provetta CryoGen® da 4 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/B	Provetta CryoGen® da 4 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BAZ	Provetta CryoGen® da 4 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BBI	Provetta CryoGen® da 4 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BG	Provetta CryoGen® da 4 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BRO	Provetta CryoGen® da 4 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BV	Provetta CryoGen® da 4 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BVI	Provetta CryoGen® da 4 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/AZ	Provetta CryoGen® da 4 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/BI	Provetta CryoGen® da 4 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/G	Provetta CryoGen® da 4 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/RO	Provetta CryoGen® da 4 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/V	Provetta CryoGen® da 4 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZ02A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL4ARIPSTS/VI	Provetta CryoGen® da 4 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEP2D	Provetta CryoGen® da 5 ml con base, tappo esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBEP2D/B	Provetta CryoGen® da 5 ml con base, tappo esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS	Provetta CryoGen® da 5 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/B	Provetta CryoGen® da 5 ml con base, tappo esterno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BAZ	Provetta CryoGen® da 5 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BBI	Provetta CryoGen® da 5 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BG	Provetta CryoGen® da 5 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BRO	Provetta CryoGen® da 5 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BV	Provetta CryoGen® da 5 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BVI	Provetta CryoGen® da 5 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/AZ	Provetta CryoGen® da 5 ml con base, tappo esterno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/BI	Provetta CryoGen® da 5 ml con base, tappo esterno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/G	Provetta CryoGen® da 5 ml con base, tappo esterno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/RO	Provetta CryoGen® da 5 ml con base, tappo esterno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/V	Provetta CryoGen® da 5 ml con base, tappo esterno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/VI	Provetta CryoGen® da 5 ml con base, tappo esterno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPS2D	Provetta CryoGen® da 5 ml con base, tappo interno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBIPS2D/B	Provetta CryoGen® da 5 ml con base, tappo interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSSN	Provetta CryoGen® da 5ml con base, tappo a vite interno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® - con bar-code a sequenza personalizzata	<i>CryoGen® Tube 5ml self standing, internal bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL5ARBIPSTS	Provetta CryoGen® da 5 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/B	Provetta CryoGen® da 5 ml con base, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BAZ	Provetta CryoGen® da 5 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BBI	Provetta CryoGen® da 5 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BG	Provetta CryoGen® da 5 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BRO	Provetta CryoGen® da 5 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BV	Provetta CryoGen® da 5 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BVI	Provetta CryoGen® da 5 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/AZ	Provetta CryoGen® da 5 ml con base, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/BI	Provetta CryoGen® da 5 ml con base, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/G	Provetta CryoGen® da 5 ml con base, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/RO	Provetta CryoGen® da 5 ml con base, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/V	Provetta CryoGen® da 5 ml con base, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/VI	Provetta CryoGen® da 5 ml con base, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS	Provetta CryoGen® da 5 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/B	Provetta CryoGen® da 5 ml, tappo interno, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BAZ	Provetta CryoGen® da 5 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BBI	Provetta CryoGen® da 5 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BG	Provetta CryoGen® da 5 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BRO	Provetta CryoGen® da 5 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BV	Provetta CryoGen® da 5 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>

Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
 | E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
 | R.N.P.P.A. IT11080P00002626





DICHIARAZIONE DI CONFORMITÀ CE / CE DECLARATION OF CONFORMITY

CL5ARIPSTS/BVI	Provetta CryoGen® da 5 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/AZ	Provetta CryoGen® da 5 ml, tappo interno AZZURRO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal LIGHT BLUE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/BI	Provetta CryoGen® da 5 ml, tappo interno BIANCO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal WHITE cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/G	Provetta CryoGen® da 5 ml, tappo interno GIALLO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal YELLOW cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/RO	Provetta CryoGen® da 5 ml, tappo interno ROSSO, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal RED cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/V	Provetta CryoGen® da 5 ml, tappo interno VERDE, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal GREEN cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/VI	Provetta CryoGen® da 5 ml, tappo interno VIOLA, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, internal VIOLET cap, graduated, with bar-code, sterile - CLEARLINE®</i>
I2D	IgenomixXera Tube 2,0 ml 2D sterili	<i>IgenomixXera Tube 2,0 ml 2D</i>



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228
| E-mail info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843
| R.N.P.P.A. IT11080P00002626