

TECHNICAL SPECIFICATIONS

Trade name	L-Shaped Spreader
Use	Disposable L-shaped sterile PS spreader for bacteriology
Material	Polystyrene PS
More	Especially handy for dispensing liquids across agar culture plates

DIMENSIONAL SPECIFICATIONS

Dimension long part L spreader	mm. 146,00
Dimension short part L spreader	mm. 32,50
Diameter	mm 5,00
General dimensional tolerance	mm. +/- 0,2



GENERAL FEATURES

CE mark	Yes
Autoclavability	No
Sterilization	Yes - γ -ray (Gamma)
Food contact suitable	Yes

CHEMICAL RESISTANCE

Categories of substances	Reference	Concen. %	T. 20°	T. 40°	T.° 60
Inorganic acid	Sulphuric acid	98	C	-	-
Organic acid	Benzoic acid	100	B	B	B
Alcohol	Ethanol	100	B	C	-
Aldehydes	Acetaldehyde	100	C	-	-
Inorganic alkali	Sodium hydroxide	50	A	A	A
Organic alkali	Aniline	100	C	-	-
Ketone	Acetone	100	-	-	-
Esther	Ethyl acetate	100	C	-	-
Hydrocarbons, halogenated	Dichloroethylene	100	C	-	-
Hydrocarbons, aromatic	Benzene	100	C	-	-
Hydrocarbons, linear	Hexane	100	C	-	-
Supersolvent	Tetrahydrofuran	100	C	-	-

Legenda (not applicable if neglected):

- **A: fair resistance;** exposure (30 days) to the chemical does not cause any damage.
- **B: sufficient resistance;** exposure causes damage of poor importance, which sometimes is only temporary.
- **C: poor resistance;** exposure to chemical is not allowed, or causes immediate permanent damage.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

PACKAGING

Case	Pcs. 1000
Pack	Single wrapped