

TECHNICAL DATA SHEET

| | |
|--------------------|----------------------|
| Part Number | 302031-100ML-HON |
| Description | Trifluoroacetic acid |

Technical specifications :

- ✓ trifluoroacetic acid chromasolv for hplc =99,0% 100ml
- ✓ CAS N°76-05-1
- ✓ Molecular weight 114.02 g/mol
- ✓ Molecular formula C₂HF₃O₂
- ✓ Density 1.48 g/cm³ (20 °C)
- ✓ Flash point 100 °C
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 100mL

- ✓ GHS pictograms GHS05 ; GHS07
- ✓ H Statements H412 ; H332 ; H314
- ✓ P Statements P260 ; P280 ; P284

Sales unit :

100mL Trifluoroacetic acid
