



PRODUCT CODE: 303117

Silver Nitrate 0.1 mol (16.987g AgNO₃) to prepare 1l of 0.1N solution

AgNO₃

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

PHYSICAL DATA: Plastic or glass ampoules containing clear colourless liquid, Miscible with water • D 1,27 • pH7 - 9 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 047-001-00-2 • VLA-ED 0,01 mg/m³



H: H314 • H411 •

P: P260 • P264 • P273 • P280 • P301+P330+P331 • P501 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P391 • P405 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

SPECIFICATIONS:

Titer at 20°C
Traceability

0,997 - 1,003
NIST

Ed.: 3 . Vig.: 03.02.2015 .

Prod.: 303117