

Honeywell

Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell Fluka
Product: 30603
CAS Numbers (All) : 1310-58-3
Molar Mass : 56.11 g/mol

Potassium hydroxide pellets, puriss. p.a., Reag. Ph. Eur.

Parameter	Specification	Units
assay	85.0-100.5	%
assay of K ₂ CO ₃	max. 1.0	%
aluminium (Al)	max. 10	ppm
calcium (Ca)	max. 0.001	%
iron (Fe)	max. 5	ppm
sodium (Na)	max. 0.5	%
nickel (Ni)	max. 0.0005	%
heavy metals (as Pb)	max. 5	ppm
chloride (Cl)	max. 100	ppm
phosphate (PO ₄)	max. 5	ppm
silicate (as SiO ₂)	max. 0.005	%
sulfate (SO ₄)	max. 100	ppm
total N	max. 0.0003	%