

Standard Analysis Certificate

PRODUCT :ISOBUTANOL RE PURO
CODE :308300
METHOD :7075

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.800 ÷ 0.804
Refractive index at 20°C	-	1.3930 ÷ 1.3990
Boiling point	°C	105.0 ÷ 109.0
Water (K.F.)	%	<=0.1
Residue on evaporation	ppm	<=50
Acidity (acetic acid)	ppm	<= 40
Assay (GLC)	%	>=99
1-Butanol	%	<= 0.4

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)