

TECHNICAL DATA SHEET

Part Number	31044-1L-HON
Description	Buffer solution pH 1.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 1,0, hydrochloric acid/potassium chloride, nist trac. and ptb 1l
- ✓ CAS N°7447-40-7
- ✓ Density 1.000 g/cm³
- ✓ Flash point no data available
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P280 ; P234

Sales unit :

1L Buffer solution pH 1.0 (20 °C)
