

TECHNICAL DATA SHEET

Part Number	31045-100ML-HON
Description	Buffer solution pH 2.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 2,0 citric acid/hcl/nacl with fungicide nist trac. and ptb 100ml
- ✓ CAS N°7647-01-0
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 100mL

- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P234 ; P280

Sales unit :

100mL Buffer solution pH 2.0 (20 °C)
