

TECHNICAL DATA SHEET

Part Number	31045-1L-HON
Description	Buffer solution pH 2.0 (20 °C)

Technical specifications :

- ✓ buffer sol. ph 2,0, citric acid/hcl/nacl, with fungicide, nist trac. and ptb 1l
- ✓ CAS N°100-56-1
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P234 ; P280

Sales unit :

1L Buffer solution pH 2.0 (20 °C)
