

TECHNICAL DATA SHEET

Part Number	31046-1L-HON
Description	Buffer solution pH 3.0 (20 °C)

Technical specifications :

- ✓ buffer sol. pH 3,0 citric acid/naoh/nacl with fungicide nist trac. and ptb 1l
- ✓ CAS N°77-92-9
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid, clear

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1L

- ✓ P Statements P280

Sales unit :

1L Buffer solution pH 3.0 (20 °C)
