

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	31048-500ML-HON
<b>Description</b>	Acetate buffer solution pH 4.6

---

## Technical specifications :

- ✓ acetate buffer sol. pH 4,6, sodium acetate / acetic acid for complexometry 500ml
- ✓ Density 1.00 g/cm<sup>3</sup>
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL
  
- ✓ P Statements P280

## Sales unit :

500mL Acetate buffer solution pH 4.6

---