

Brand: Honeywell Fluka  
Product: 31232

Iron(III) chloride hexahydrate puriss. p.a., Reag. Ph. Eur.

Parameter	Specification	Units
assay	min. 99	%
free acid (as HCl)	max. 0.2	%
arsenic (As)	max. 0.0005	%
copper (Cu)	max. 0.002	%
iron (Fe 2+)	max. 0.005	%
lead (Pb)	max. 0.002	%
zinc (Zn)	max. 0.002	%
sub.not precip.b.NH4OH(SO4)	max. 0.05	%
free chlorine (Cl)	max. 0.001	%
phosphate (PO4)	max. 0.002	%
sulfate (SO4)	max. 0.005	%
total N	max. 0.001	%