

Brand: Honeywell Fluka
Product: 31263

Potassium nitrate puriss. p.a., Reag. ACS, Reag. ISO, Reag. Ph. Eur.

Parameter	Specification	Units
assay	min. 99.0	%
pH (5 %, 25 °C)	5.0-8.0	
water insoluble matter	max. 0.005	%
ammonium (NH ₄)	max. 0.001	%
calcium (Ca)	max. 0.001	%
copper (Cu)	max. 0.0005	%
iron (Fe)	max. 2	ppm
magnesium (Mg)	max. 0.002	%
sodium (Na)	max. 0.005	%
lead (Pb)	max. 0.0005	%
zinc (Zn)	max. 0.0005	%
heavy metals (as Pb)	max. 5	ppm
chloride (Cl)	max. 0.001	%
iodate (IO ₃)	max. 5	ppm
nitrite (NO ₂)	max. 0.001	%
phosphate (PO ₄)	max. 5	ppm
sulfate (SO ₄)	max. 0.002	%