

TECHNICAL DATA SHEET

Part Number	319252-2L-HON
Description	Phosphate buffer solution

Technical specifications :

- ✓ phosphate buffer solution, ph 7,2 2l
- ✓ CAS N°1310-73-2
- ✓ Density 0.995 g/cm³ (25 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 2L

- ✓ P Statements P280

Sales unit :

2L Phosphate buffer solution
