

TECHNICAL DATA SHEET

Part Number	319252-500ML-HON
Description	Phosphate buffer solution

Technical specifications :

- ✓ phosphate buffer solution, ph 7,2 500ml
- ✓ CAS N°7778-77-0
- ✓ Density 0.995 g/cm³ (25 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL

- ✓ P Statements P280

Sales unit :

500mL Phosphate buffer solution
