

## Stomacher® Bags - Strainer Bags

### Product Codes by Sampling Volume

5-80ml: **BA6040/STR, BA6040/STR/CLR, BA6040/STR/DBL**  
80-400ml: **BA6041/STR, BA6041/LSTR, BA6141/STR**


### TECHNICAL SPECIFICATIONS

- Stomacher® blender bag, with addition of strainer element
- Outer bag – Virgin grade, Polyolefin complex. Flexible and transparent
- Strainer bags – Virgin grade, microfiltered Polyolefin complex. Flexible and opaque
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenisation of the sample
- Compatible with any paddle lab blender
- Suitable for food contact. Not designed to be used for food packaging
- Irradiated sterile in sealed sachets, with batch traceable certificate
- Storage: in a temperate and dry place, away from humidity and light
- Available in 80ml and 400ml

### Stomacher® 80 Biomaster Strainer Bag


P/N. **BA6040/STR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)




P/N. **BA6040/STR/CLR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)
- Additional closure tab, consisting of wire reinforced tape



P/N. **BA6040/STR/DBL**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)



### Seward Ltd

Technology Centre  
Easting Close  
Worthing, BN14 8HQ  
United Kingdom  
info@seward.co.uk

### Seward Lab Systems Inc

155 Keyland Court  
Bohemia  
NY 11716  
USA  
info@sewardusa.com


### FermionX GmbH

Münstertäler Str. 12,  
Eschbach  
DE-79427  
Germany  
info@seward.co.uk

## Stomacher® 400 Strainer Bags


### P/N. BA6041/STR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Straight bottom seal
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)




### P/N. BA6041/LSTR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Lateral side strainer arrangement
- Interwoven microfiber filter
- Straight bottom seal
- Pack size: 500 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



### P/N. BA6141/STR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Full strainer bag arrangement
- Microfiltered, 0.5-0.6mm pore size
- Round bottom seal
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



Stomacher® is registered trademark of Seward Ltd.

SEWARD is a part of the FermionX Group. FermionX Ltd reserves the right to change and improve specifications of its products without notice. Please contact [info@seward.co.uk](mailto:info@seward.co.uk) or visit [www.seward.co.uk](http://www.seward.co.uk) for more information.

#### Seward Ltd

Technology Centre  
Easting Close  
Worthing, BN14 8HQ  
United Kingdom  
[info@seward.co.uk](mailto:info@seward.co.uk)

#### Seward Lab Systems Inc

155 Keyland Court  
Bohemia  
NY 11716  
USA  
[info@sewardusa.com](mailto:info@sewardusa.com)

#### FermionX GmbH

Münstertäler Str. 12,  
Eschbach  
DE-79427  
Germany  
[info@seward.co.uk](mailto:info@seward.co.uk)