

Stomacher® Bags - Standard

Product Codes by Sampling Volume

250ul to 80ml: **BA6040**

80-280ml: **BA6241, BA6341 (Eco bags)**

80-400ml: **BA6041, BA6141,**

400-4500ml: **BA6042, BA6042/HD, BA6042/ACL**


TECHNICAL SPECIFICATIONS

- Stomacher® blender bag, no strainer or closure addition
- Virgin grade, Polyolefin complex
- Flexible and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenisation of the sample
- Compatible with any paddle lab blender
- Suitable for food contact. Not designed to be used for food packaging
- Irradiated sterile in sealed sachets, with batch traceable certificate
- Storage: in a temperate and dry place, away from humidity and light
- Available in 80ml, 400ml, 3500 mL

Stomacher® 80 Biomaster Standard Bag

P/N. **BA6040**


- Max Blending Volume: 250ul to 80ml
- Recommended blending Volume: 80ml
- Pack size: 2000 Bags (20 x 50 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)



Stomacher® 400 Eco Bag (Path.)

P/N. **BA6241**


- Max Blending Volume: 250ml
- Recommended blending Volume: 250ml
- Straight bottom seal
- Pack size: 500 Bags (10 x 50 sachets)
- Bag Dimensions: 177x200ml (65µm thickness)



Stomacher® 400 Eco Bag (TVC)

P/N. **BA6341**

- Max Blending Volume: 200ml
- Recommended blending Volume: 100ml
- Straight bottom seal
- Pack size: 500 Bags (10 x 50 sachets)
- Bag Dimensions: 155x200ml (65µm thickness)



Seward Ltd

Technology Centre
Easting Close
Worthing, BN14 8HQ
United Kingdom
info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
Bohemia
NY 11716
USA
info@sewardusa.com

FermionX GmbH

Münstertäler Str. 12,
Eschbach
DE-79427
Germany
info@seward.co.uk

Stomacher® 400 Classic Standard Bag

P/N. **BA6041**

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Straight bottom seal
- Pack size (BA6041): 500 Bags (10 x 50 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)

Stomacher® 400 Circulator Standard Bag

P/N. **BA6141**

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Curved bottom seal
- Pack size (BA6141): 500 Bags (10 x 50 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)

Stomacher® 3500 Jumbo Standard Bag

P/N. **BA6042**

- Max Blending Volume: 4500ml
- Recommended blending Volume: 3750ml
- Straight bottom seal
- Pack size (BA6042): 250 Bags (5 x 50 sachets)
- Bag Dimensions: 380x510ml (65µm thickness)

Stomacher® 3500 Jumbo Heavy-Duty Bag





P/N. **BA6042/HD**

- Max Blending Volume: 4500ml
- Recommended blending Volume: 3750ml
- Straight bottom seal
- Pack size (BA6042/HD): 250 Bags (5 x 50 sachets)
- Bag Dimensions: 380x510ml (100µm thickness)

Stomacher® 3500 Autoclavable Blender Bag

P/N. **BA6042/ACL**

- Max Blending Volume: 4500ml
- Recommended blending Volume: 3750ml
- Straight bottom seal
- Pack size (BA6042/ACL): 500 Bags (10 x 50 sachets)
- Bag Dimensions: 380x510ml (100µm thickness)



Stomacher® is registered trademark of Seward Ltd.

SEWARD is a part of the FermionX Group. FermionX Ltd reserves the right to change and improve specifications of its products without notice. Please contact info@seward.co.uk or visit www.seward.co.uk for more information.

Seward Ltd

Technology Centre
Easting Close
Worthing, BN14 8HQ
United Kingdom
info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
Bohemia
NY 11716
USA
info@sewardusa.com

FermionX GmbH

Münstertäler Str. 12,
Eschbach
DE-79427
Germany
info@seward.co.uk