

## TECHNICAL DATA SHEET

---

<b>Part Number</b>	3193178
<b>Description</b>	Bag rack



---

### Benefits :

- ✓ Seward bag racks provide the perfect storage for Stomacher® bag sample handling, helping to provide a clean and productive working environment.
- ✓ Perfect for loading closure bags preenrichments into the incubator
- ✓ Made from high grade plastic polymer
- ✓ Easily dismantled and cleaned with most lab disinfectants (not autoclavable)

### Technical specifications :

- ✓ Dimensions : 8 x 10.5 x 19 cm
- ✓ Weight : 0.4 kg

### Sales unit :

1 Bag rack

---