

*Flexible bench or
under-shelf storage*

Flip & Grip™ Pipette Holder

The Flip & Grip™ pipette holder provides flexibility for changing storage needs. When used upright, the stand provides stable bench-top storage for up to 12 single-channel pipettes.

If space is tight or overhead storage is preferred, just unscrew the top knob and reverse the holding ring.

The unit can now be mounted via a screw in the base to any appropriate surface (such as a shelf or undercounter space) to provide hanging storage.

Rubber coated grips hold each pipette firmly in place without scratching, and with enough distance between units to ensure no contact between pipette tips. A weighted base provides stability during carousel rotation, even when unit is unevenly loaded.

The stand holds up to 12 single-channel pipettes or up to 3 single/3 multichannel pipettes (six spaces for pipette barrels and six hooks). It will accommodate most single-channel pipette brands.

Bench-top (upright) dimensions:
332 mm H x 193 mm W.

Length when shelf-mounted (upside-down): 187 mm.




Mount base under shelf or counter



**Rack rotates
for easy pipette
selection**

Flexible bench or
under-shelf storage

Flip & Grip™ Pipette Holder



*Pipettes can
be held in the
rubber grips
and/or the
plunger hooks
available*

Available in
3 colors

part no.	color
HS120026	Blue
HS120056	Lime
HS120057	Magenta













Heathrow Scientific offers a full line of Pearl™ pipettes
Precision pipetting at your fingertips



Pearl™ pipettes were specifically engineered for ease of use and comfort. Our engineers reviewed and re-designed each functional part to create the simplest, most comfortable option available on the market.

The Pearl™ family of air-displacement pipettes includes ten adjustable volume instruments that cover a complete range of volumes from 0.1 µL to 10,000 µL.

color	part no.	volume range
	HS120097	0.1–2µL
	HS120098	0.5–10µL
	HS120099	2–20µL
	HS120100	5–50µL
	HS120101	10–100µL
	HS120102	20–200µL
	HS120103	30–300µL
	HS120104	100–1000µL
	HS120105	500–5000µL
	HS120106	1000–10000µL

Heathrow Scientific, LLC
620 Lakeview Parkway, Vernon Hills IL 60061 USA
Tel.: +1-800-741-4597 | Fax: +1-847-816-5072 | E-mail: info@heatsci.com