

## TECHNICAL DATA SHEET

---

<b>Part Nbr</b>	320126
<b>Designation</b>	81-well microtube storage box



---

### Description :

- ✓ Made from polypropylene
- ✓ Easy open, friction fit lid
- ✓ Designed for quick visual identification with uniform wells and location identification
- ✓ Stackable, saves valuable space in the lab
- ✓ Store down to -80°C - Must use tubes capable of being stored at the required temperature of the user
- ✓ Autoclavable
- ✓ Tubes: 1,5/2,0 ml
- ✓ Wells: 81 – wells Ø: 12,3 mm
- ✓ Color : natural
- ✓ Dimensions: 13 x 13 x 4,7 cm

### Sales unit :

5x 81-well microtube storage box

---