

TECHNICAL DATA SHEET

Part Number	32043-1EA-HON
Description	Sulfuric acid concentrate

Technical specifications :

- ✓ sulfuric acid concentrate, for 1l standard solution, 0,05 m h₂so₄ (0,1n)
- ✓ CAS N°7664-93-9
- ✓ Molecular weight 98.08 g/mol
- ✓ Molecular formula H₂O₄S
- ✓ Density 1.05 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05
- ✓ H Statements H319 ; H315 ; H290
- ✓ P Statements P280 ; P234

Sales unit :

1 Each Sulfuric acid concentrate
