



# Premium Plus Slide Boxes

Modern design for daily use

**NEW!**





- a. Slide rail design ensure a secure fit allowing easy placement and removal of slides
- b. Centrally positioned location numbering grid enables simple slide organization with easy viewing of location numbers
- c. Entirely made from ABS and Polypropylene makes the boxes robust and eliminates potential rusting of metal parts

- d. Sliding lock provides secure storage with easily identified lock and unlock symbols. 2 mm zip tie can be used for more secure storage.
- e. Inventory card is removable and replaceable for trouble-free writing on flat surfaces. Inventory cards can be purchased separately in pks/10
- f. Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement.

- g. EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- h. Multiple labeling areas, a large area on the lid of the box with "project" and "date" molded in and the back of the index card can be used as the lid is transparent

- i. Vibrant translucent colors are easy to identify and enable coding options, assorted 5/pk available
- j. Stylish wave surface texture is easy to hold, increasing safety of sample transport
- k. Stackable, molded feet allow secure horizontal stacking of multiple units for safe transport or storing

- l. Flexible storing, the units can be stored vertically or horizontally helping to efficiently use space available
- m. Unique side storage can accommodate a pen, desiccant or other small items needed for working
- n. Holds 100 - 26 x 75 mm slides

Patent Pending

DOMINIQUE SCHLIER SAS

# NEW!



## 100-Place Premium Plus Slide Box

| Item No. | Color     | Lining   | Hinge         | L x W x H<br>in | L x W x H<br>cm   | UOM                                   |
|----------|-----------|----------|---------------|-----------------|-------------------|---------------------------------------|
| 120576   | Blue      | EVA Foam | Sliding Latch | 8.2 x 7.7 x 1.3 | 20.7 x 19.5 x 3.2 | 5/pk                                  |
| 120577   | Orange    |          |               |                 |                   | 4/pks per case (shrink wrapped packs) |
| 120578   | Green     |          |               |                 |                   |                                       |
| 120579   | Red       |          |               |                 |                   |                                       |
| 120592   | Purple    |          |               |                 |                   |                                       |
| 120593   | Assorted* |          |               |                 |                   |                                       |

\*Assorted colors (Purple, Blue, Green, Red, Orange)

## Inventory Card

| Item No. | L x H<br>in | L x H<br>cm | UOM                                      |
|----------|-------------|-------------|--|
| 120618   | 8.5 x 11    | 21.6 x 28   | 10/PK (Polybag with Label)50 PK per case |