



# Premium Plus Slide Boxes

Modern design for daily use **NEW!**





- a. Slide rail design ensure a secure fit allowing easy placement and removal of slides
- b. Centrally positioned location numbering grid enables simple slide organization with easy viewing of location numbers
- c. Entirely made from ABS and Polypropylene makes the boxes robust and eliminates potential rusting of metal parts

- d. Sliding lock provides secure storage with easily identified lock and unlock symbols. 2 mm zip tie can be used for more secure storage.
- e. Inventory card is removable and replaceable for trouble-free writing on flat surfaces. Inventory cards can be purchased separately in pks/10
- f. Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement.

- g. EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- h. Multiple labeling areas, a large area on the lid of the box with "project" and "date" molded in and the back of the index card can be used as the lid is transparent

- i. Vibrant translucent colors are easy to identify and enable coding options, assorted 5/pk available
- j. Stylish wave surface texture is easy to hold, increasing safety of sample transport
- k. Stackable, molded feet allow secure horizontal stacking of multiple units for safe transport or storing

- l. Flexible storing, the units can be stored vertically or horizontally helping to efficiently use space available
- m. Unique side storage can accommodate a pen, desiccant or other small items needed for working
- n. Holds 100 - 26 x 75 mm slides

Patent Pending

DOMINIQUE SCHLIER SAS



# NEW!



## 100-Place Premium Plus Slide Box

Item No.	Color	Lining	Hinge	L x W x H in	L x W x H cm	UOM
120576	Blue	EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2	5/pk
120577	Orange					4/pks per case (shrink wrapped packs)
120578	Green					
120579	Red					
120592	Purple					
120593	Assorted*					

\*Assorted colors (Purple, Blue, Green, Red, Orange)

## Inventory Card

Item No.	L x H in	L x H cm	UOM
120618	8.5 x 11	21.6 x 28	10/PK (Polybag with Label)50 PK per case