

# Mini LED Light Box

The small and slim footprint saves space on the bench-top and can be conveniently stored in a drawer.

- Ideal for quick screening of 96-well plates, Petri dishes, electrophoresis gels, slides, and transparencies
- LED light source stays cool and has longer life (10,000 hours) and low power consumption
- Viewing area 125 x 100 mm (5" x 4")
- Powered by six 1.5V AAA

## Mini Magnifier

The free standing magnifier enlarges objects for enhanced viewing

- Molded lens swivels to adjust viewing angle for optimal view
- 1.75x magnification

Item No.	Description	L x W x H in	L x W x H cm	UOM
120661	Mini LED Light Box*	6.2 x 5.9 x 0.6	15.8 x 15 x 1.4	1/ea
HS15092D	Mini Magnifier	6.3 x 5.3 x 8.1	16.1 x 13.5 x 20.7	1/ea, 24ea/cs

*\*Batteries included*



Both the Mini LED Light Box and Mini Magnifier are the perfect accessories with the eCount Colony Counter!