

# Reagent Reservoirs

Reagent Reservoirs by Heathrow Scientific offer a fresh and modern design with enhanced features to streamline pipetting procedures.

- A full line of sizes are available in both sterile and non-sterile
- Made from Clear PVC or White Polystyrene
- Trough design with a “v” shaped bottom helps ensure that all the reagent is accessible and can be easily aspirated
- Thick wall construction and wide base adds stability and durability making the reservoir easy to handle without product deformation which results in less spillage
- Designed to stay flat on the benchtop and not warp
- Sterile versions are sealed in bags and electron beam sterilized. Dnase, Rnase and pyrogen free
- Sealed in easy-to-open bags. No scissors needed to open
- Quality certificate available on request
- Disposable
- Design Patent Pending

## Stackable

Reservoirs can be conveniently stacked with a unique nesting design that enables easy separation.




## Divided Reservoir Design

allows pipetting of two different reagents or maximizes recovery if using only single channel and small amount of reagent.



### Graduation Markings

on the top edge and inner walls allows measurements of liquid to be easily identified and reduces wastage of expensive reagents.



### Pour-off Spouts

Improved design of pour-off spouts on all four corners enables easy spill-free pouring and reduces chance of dripping.



120647, 120648  
120649



120650



120643, 120644  
120645



120646, 120690



120637



120638, 120639



120651, 120652



120640, 120641  
120642

Item No.	Volume mL	Material	Sterility	Interior Dimensions L x W x H in	Interior Dimensions L x W x H cm	UOM
120637	5 mL	Polystyrene	Yes	2.8 x 0.8 x 0.7	7.2 x 2.1 x 1.8	1 ea/bag, 50 bags/pk, 10 pks/cs
120638	10 mL	Polystyrene	No	2.8 x 0.8 x 0.7	7.2 x 2.1 x 1.8	100 ea/pk, 6 pks/cs
120639			Yes			1 ea/bag, 100 bags/pk, 6 pks/cs
120640	25 mL	Polystyrene	Yes	4 x 1 x 0.8	10.2 x 2.5 x 2.1	5 ea/bag, 40 bags/pk, 10 pks/cs
120641			Yes			1 ea/bag, 100 bags/pk, 6 pks/cs
120642			No			100 ea/pk, 6 pks/cs
120643			Yes			5 ea/bag, 40 bags/pk, 10 pks/cs
120644	50 mL	Polystyrene	Yes	4 x 1.1 x 1.2	10.2 x 2.8 x 3	1 ea/bag, 100 bags/pk, 6 pks/cs
120645			No			100 ea/pk, 6 pks/cs
120646			Yes			5 ea/bag, 20 bags/pk, 6 pks/cs
120690			No			100 ea/pk, 6 pks/cs
120647	100 mL	Polystyrene	Yes	4.3 x 1.9 x 1.3	10.9 x 4.8 x 3.2	5 ea/bag, 40 bags/pk, 6 pks/cs
120648			Yes			1 ea/bag, 100 bags/pk, 6 pks/cs
120649			No			100 ea/pk, 6 pks/cs
120650			No			100 ea/pk, 6 pks/cs
120651	25 mL (10 and 15 mL divided)	Polystyrene	Yes	4 x 1 x 0.8	10.2 x 2.5 x 2.1	5 ea/bag, 40 bags/pk, 6 pks/cs
120652			Yes			1 ea/bag, 100 bags/pk, 6 pks/cs