

Product Highlights

- Accommodates all brands of 13 and 16 mm round, flat and conical bottom blood collection & culture tubes
- Six way connection, link different or the same well sizes, end to end or side by side, in any orientation
- Customizable racks can be linked together, or used as a standalone
- Assorted color pack available for color coding options
- Racks can be used in water baths and placed in freezers

Holds
Multiple Tube Sizes
without multiple
racks!



Assorted Color Pack
available for color coding options

Blood Tube Brands (include & not limited to):


- Becton Dickinson
- Greiner
- Covidien
- Thermo Scientific™
- Applied Biosystems
- Sarstedt

Recommended Tube Dimensions:

- 13 x 75 mm
- 13 x 100 mm
- 16 x 100 mm
- 16 x 125 mm
- 17 x 120 mm


Tube Volume Guide:

- 3.5 mL
- 7.0 mL
- 10 mL
- 15 mL



Customizable

Racks can slide together to create an infinite number of combinations




Clinical 4-Way Tube Rack

- Accommodates all brands of 13 and 16 mm round, flat and conical bottom blood collection & culture tubes
- Only 4-way rack available with a six way connection, link different or the same well sizes, end to end or side by side, in any orientation
- Customize the rack to fit the number of wells required, infinite units can be linked together or use one standalone unit
- Uniformly spaced wells help eliminate pipetting errors which enables more accurate results
- Molded alphanumeric identification for quick location and retrieval of tubes
- Robust, made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezers
- Assorted color pack available for color coding options
- Patent Pending
- **Autoclavable**



Molded Alphanumeric Identification
for quick location and retrieval of tubes



Tube Bottoms Nest
in well pockets keeping them from going thru the rack and off the bench top




Versatile

Accommodates tubes with round, flat and conical bottoms

↑
Holds 15 mL
Conical Tubes!

Item No.	Color	Item No.	Color
120684	Blue	120687	Natural
120685	Purple	120683	Assorted*
120686	Green		

*Assorted colors (Blue, Purple, Green, Red, Natural)



Individual Rack Specifications

Side	Wells	Well Dia. mm	L x W x H in	cm	UOM
A	12	13	9.1 x 2.8 x 3.8	23.2 x 7.0 x 9.5	5 ea/pk, 6 pk/cs
B	18	13			
C	12	17			
D	18	17			