


## Product Highlights


- Chemically-resistant HDPE Coated Wire Racks
- Racks can be used between temperatures of -40°C and 121°C
- Holds all brands and types of conical, round or flat bottom, glass or plastic, tubes and/or bottles
- Square wells and three deck levels help keep tubes secure
- Autoclavable



\*HS120089  
6x12, 72 Wells

### Fits 10-13mm Tubes

| Item No.  | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm      | UOM            |
|-----------|-----------------|------------|-------|---------------|-----------------|-------------------|----------------|
| HS120089* | Blue            | 6x12       | 72    | 13            | 8.1 x 4.1 x 2.6 | 20.5 x 10.5 x 6.5 | 1 ea, 48 ea/cs |
| 120759    | <b>NEW</b> Blue | 9x12       | 108   | 13            | 8 x 6.4 x 2.6   | 20.4 x 16.1 x 6.5 | 1 ea, 36 ea/cs |



\*120761  
4x12, 48 Wells

### Fits 13-16mm Tubes

| Item No. | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm      | UOM            |
|----------|-----------------|------------|-------|---------------|-----------------|-------------------|----------------|
| 120760   | <b>NEW</b> Blue | 4x10       | 40    | 16            | 7.9 x 3.3 x 2.4 | 20 x 8.3 x 5.9    | 1 ea, 56 ea/cs |
| 120761*  | <b>NEW</b> Blue | 4x12       | 48    | 16            | 9.4 x 3.3 x 2.4 | 23.7 x 12.7 x 6.5 | 1 ea, 56 ea/cs |
| HS23072  | Blue            | 6x12       | 72    | 16            | 9.5 x 5 x 2.6   | 24.1 x 18.8 x 6.5 | 1 ea, 24 ea/cs |
| 120763   | <b>NEW</b> Blue | 9x12       | 108   | 16            | 9.5 x 7.4 x 2.6 | 24.1 x 18.8 x 6.5 | 1 ea, 24 ea/cs |

### Fits 16-18mm Tubes – 15mL


| Item No. | Color | Tube Array | Wells | Well Diam. mm | L x W x H in | L x W x H cm     | UOM            |
|----------|-------|------------|-------|---------------|--------------|------------------|----------------|
| HS23215* | Blue  | 3x5        | 15    | 18            | 5 x 3 x 3    | 12.7 x 7.6 x 7.6 | 1 ea, 96 ea/cs |



\*HS23215  
3x5,  
15 Wells

### Fits 16-20mm Tubes – 15mL

| Item No. | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm | UOM            |
|----------|-----------------|------------|-------|---------------|-----------------|--------------|----------------|
| 120764   | <b>NEW</b> Blue | 4x10       | 40    | 20            | 9.5 x 3.9 x 3.1 | 24 x 9 x 8   | 1 ea, 48 ea/cs |
| 120765   | <b>NEW</b> Blue | 8x10       | 80    | 20            | 9.5 x 7.6 x 3.1 | 24 x 19 x 8  | 1 ea, 24 ea/cs |



\*120768  
8x10, 80 Wells

### Fits 20-24mm Tubes

| Item No. | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in     | L x W x H cm  | UOM            |
|----------|-----------------|------------|-------|---------------|------------------|---------------|----------------|
| 120766   | <b>NEW</b> Blue | 4x10       | 40    | 24            | 10.2 x 4.2 x 3.2 | 26 x 10.7 x 8 | 1 ea, 24 ea/cs |

### Fits 20-25mm Tubes

| Item No. | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in     | L x W x H cm    | UOM            |
|----------|-----------------|------------|-------|---------------|------------------|-----------------|----------------|
| 120769   | <b>NEW</b> Blue | 3x8        | 24    | 25            | 10.7 x 4.1 x 3.2 | 27.3 x 10.5 x 8 | 1 ea, 24 ea/cs |
| 120767   | <b>NEW</b> Blue | 4x10       | 40    | 25            | 11.4 x 4.7 x 3.1 | 29 x 11.8 x 8   | 1 ea, 24 ea/cs |
| 120768*  | <b>NEW</b> Blue | 8x10       | 80    | 25            | 11.4 x 9.1 x 3.1 | 29 x 23 x 8     | 1 ea, 12 ea/cs |




\*HS232461  
2x4, 8 Wells

## Fits 25-30mm Tubes - 50mL


| Item No.  | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in     | L x W x H cm      | UOM            |
|-----------|-----------------|------------|-------|---------------|------------------|-------------------|----------------|
| HS232461* | Blue            | 2x4        | 8     | 30            | 5.8 x 3.2 x 3.5  | 14.7 x 8 x 8.9    | 1 ea, 96 ea/cs |
| 120770    | <b>NEW</b> Blue | 6x8        | 48    | 30            | 10.9 x 8.3 x 3.9 | 27.7 x 20.9 x 9.8 | 1 ea, 8 ea/cs  |

## Fits 30-40mm Tubes

| Item No. | Color           | Tube Array | Wells | Well Diam. mm | L x W x H in     | L x W x H cm      | UOM            |
|----------|-----------------|------------|-------|---------------|------------------|-------------------|----------------|
| 120771*  | <b>NEW</b> Blue | 3x8        | 24    | 40            | 14.1 x 5.5 x 3.9 | 35.7 x 13.9 x 9.8 | 1 ea, 12 ea/cs |
| 120772   | <b>NEW</b> Blue | 6x8        | 48    | 40            | 14 x 10.6 x 3.8  | 35 x 26 x 9.8     | 1 ea, 8 ea/cs  |



\*120771  
3x8, 24 Wells




\*HS120086  
5x10, 50 Wells

## Fits 23-25mm Bottles - Bijou Bottles


| Item No.  | Color | Tube Array | Wells | Well Diam. mm | L x W x H in     | L x W x H cm      | UOM            |
|-----------|-------|------------|-------|---------------|------------------|-------------------|----------------|
| HS120085  | Blue  | 2x6        | 12    | 25            | 6.4 x 2.5 x 1.3  | 16.4 x 6.4 x 3.2  | 1 ea, 96 ea/cs |
| HS120086* | Blue  | 5x10       | 50    | 25            | 10.7 x 5.7 x 1.3 | 27.1 x 14.5 x 3.2 | 1 ea, 48 ea/cs |

## Fits 30-33mm Bottles - Universal & McCarthy Bottles

| Item No.  | Color | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm    | UOM            |
|-----------|-------|------------|-------|---------------|-----------------|-----------------|----------------|
| HS120087* | Blue  | 2x6        | 12    | 33            | 8.3 x 3.1 x 2.4 | 21.2 x 8 x 6    | 1 ea, 72 ea/cs |
| HS120088  | Blue  | 5x10       | 50    | 33            | 14 x 7.4 x 2.4  | 35.5 x 18.9 x 6 | 1 ea, 24 ea/cs |



\*HS120087  
2x6, 12 Wells



\*HS120091  
4x6, 24 Wells

## Fits 55-58mm Bottles - Urine Containers

| Item No.  | Color | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm    | UOM            |
|-----------|-------|------------|-------|---------------|-----------------|-----------------|----------------|
| HS120091* | Blue  | 4x6        | 24    | 58            | 15.5 x 10.3 x 2 | 39.3 x 26.2 x 5 | 1 ea, 12 ea/cs |










## Fits 55-60mm Bottles - Centrifuge Bottles





| Item No. | Color | Tube Array | Wells | Well Diam. mm | L x W x H in    | L x W x H cm       | UOM            |
|----------|-------|------------|-------|---------------|-----------------|--------------------|----------------|
| HS23224* | Blue  | 2x3        | 6     | 60            | 8.1 x 5.3 x 4.1 | 20.5 x 13.5 x 10.4 | 1 ea, 24 ea/cs |
















\*HS23224  
2x3, 6 Wells

TUBE STYLE SELECTOR

|                               |   |   |   |   |   |   |   |   |   |
|-------------------------------|---|---|---|---|---|---|---|---|---|
| Example Tube Style            |  |  |  |  |  |  |  |  |  |
| Example Tube/Bottle/Vial Type | Blood Tubes<br>Culture Tubes<br>Test Tubes  | Blood Tubes<br>Culture Tubes<br>Test Tubes  | Test Tubes<br>Culture Tubes<br>15mL Tube  | Test Tubes<br>Culture Tubes<br>15mL Tube  | Test Tubes<br>Culture Tubes   | Test Tubes<br>Culture Tubes   | Test Tubes<br>Culture Tubes<br>50mL Tube  | Test Tubes<br>Culture Tubes   |   |
| Well Tolerance                | 10-13 mm  | 13-16mm   | 16-18mm   | 16-20mm   | 20-24mm   | 20-25mm   | 25-30mm   | 30-40mm   |   |
| Rack Codes                    | HS120089<br>120759  | 120761<br>HS23072<br>120763<br>120760   | HS23215   | 120764<br>120765  | 120766  | 120767<br>120768<br>120769  | 120770<br>HS232461  | 120771<br>120772  |   |

|                               |   |   |  |   |
|-------------------------------|---|---|--|---|
| Example Tube Style            |  |          |  |  |
| Example Tube/Bottle/Vial Type | Bijou Bottle<br>Scintillation Vial  | McCarthy Bottle<br>Universal Bottle<br>Scintillation Vial<br>Culture Bottle<br>Serum Vial | Urine Containers<br>Centrifuge Bottles<br>Culture Bottles<br>Storage Bottles       | Centrifuge Bottles<br>Culture Bottles<br>Storage Bottles                            |
| Well Tolerance                | 23-25mm   | 30-33mm   | 55-58mm  | 55-60mm   |
| Rack Codes                    | HS120085<br>HS120086  | HS120087<br>HS120088  | HS120091   | HS23224   |

TUBE ARRAY EXAMPLES

|   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|
| 2x3   | 2x4   | 2x6   | 3x5   | 3x8   | 4x6   | 4x10  |
|  |  |  |  |   |  |  |
| 4x12  | 5x10  | 6x8   | 6x12  | 8x10  | 9x12  |   |
|  |  |  |  |  |  |   |