

Honeywell Specialty Chemicals Seelze GmbH  
Wunstorferstrasse 40  
Seelze, 30926  
Telefon: +49 5137 999-0  
[lab.honeywell.com](http://lab.honeywell.com)

Brand: Honeywell Fluka  
Product: 32301

Molar Mass : 77.08 g/mol  
CAS Numbers (All) : 631-61-8

Ammonium acetate puriss. p.a., Reag. ACS, Reag. Ph. Eur.

Parameter	Specification	Units
assay	min. 98	%
water insoluble matter	max. 0.005	%
water (Karl Fischer)	max. 2	%
sulfated ash	max. 0.01	%
pH (5 %, 25 °C)	6.7-7.3	
calcium (Ca)	max. 0.001	%
cadmium (Cd)	max. 0.0002	%
cobalt (Co)	max. 0.0005	%
copper (Cu)	max. 0.0002	%
chromium (Cr)	max. 0.0005	%
iron (Fe)	max. 0.0005	%
potassium (K)	max. 0.005	%
magnesium (Mg)	max. 0.0002	%
manganese (Mn)	max. 0.0005	%
sodium (Na)	max. 0.005	%
nickel (Ni)	max. 0.0005	%
lead (Pb)	max. 0.0002	%
zinc (Zn)	max. 0.0002	%
heavy metals (as Pb)	max. 0.0005	%
chloride (Cl)	max. 0.0005	%
nitrate (NO <sub>3</sub> )	max. 0.001	%
sulfate (SO <sub>4</sub> )	max. 0.001	%
KMNO <sub>4</sub> red. matter(as HCOOH)	max. 0.005	%