

Tray Plate (#30001)

1. Description

SPL's Tray Plates provide a larger surface area and easy handling compared to conventional round type dishes. The external dimensions of Tray Plates are identical to standard SPL Cell Culture Plates for broader applications

2. Features

- ❖ Designed to prevent cross contamination
- ❖ Effective gas exchange lid inner design

3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase free	:	Passed
DNA free	:	Passed

4. Product Range

Cat.No.	Material	External Dimension w x l x h (mm)	Well Dimension w x l x h (mm)	Growth Area (cm ²)	Surface Treatment	Sterile	Packaging
30001	PS	127.95 x 85.5 x 16.2	105.47 x 69.45 x 11.6	80	+	+	1/50

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photograph




7. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging
Lid	Plate		
Length: 84.5 mm Width: 126.5 mm Height: 9.00 mm	Length: 85.5 mm Width: 127.95 mm Height: 16.20 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> • Pieces/Sleeve: 10 • Pieces/Case: 50

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-TP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, November 20th, 2018