

## Aspiration Pipette

(#94001/94002/94005/94010)

### 1. Description

SPL Aspiration Pipettes are useful for quick collection and discardment of liquid samples.

### 2. Features

- ❖ No graduation scales
- ❖ No plugs

### 3. Functional Test

Sterility	:	Passed
Cytotoxicity	:	Passed
Non-Pyrogenic	:	Passed
RNase/DNase free	:	Passed
DNA free	:	Passed

### 4. Product Range

Cat.No.	Material	Working Volume (ml)	Sterile	Packaging
94001	PS	1.00	+	1 / 200 / 800
94002	PS	2.00	+	1 / 150 / 600
94005	PS	5.00	+	1 / 100 / 400
94010	PS	10.00	+	1 / 100 / 400

### 5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

### 6. Product Photographs



SPL-TDS-ASP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
© SPL Life Sciences Co., Ltd., Version 2, July 19<sup>th</sup>, 2016

## 7. Specifications &amp; Technical Drawing

(1) 1 ml (# 94001)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 277.00 mm Diameter: 4.60 mm Volume: 1.00 ml	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 1</li> <li>• Pieces/Inbox: 200</li> <li>• Pieces/Case: 800</li> </ul>

- Technical Drawing



(2) 2 ml (# 94002)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 277.00 mm Diameter: 6.40 mm Volume: 2.00 ml	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 1</li> <li>• Pieces/Inbox: 150</li> <li>• Pieces/Case: 600</li> </ul>

- Technical Drawing




(3) 5 ml (# 94005)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 338.00 mm Diameter: 8.30 mm Volume: 5.00 ml	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 1</li> <li>• Pieces/Inbox: 100</li> <li>• Pieces/Case: 400</li> </ul>

- Technical Drawing




(4) 10 ml (# 94010)

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 338.00 mm Diameter: 9.90 mm Volume: 10.00 ml	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 1</li> <li>• Pieces/Inbox: 100</li> <li>• Pieces/Case: 400</li> </ul>

- Technical Drawing

[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&amp;D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-ASP

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
© SPL Life Sciences Co., Ltd., Version 2, July 19<sup>th</sup>, 2016