

## Microtube Rack

(#61080)

### 1. Description

SPL Microtube Racks are designed for handling and storage of Microcentrifuge Tubes.

### 2. Features

- ❖ Numeric labeling
- ❖ 5 colors: White, Blue, Green, Red, Purple
- ❖ Autoclavable

### 3. Product Range

Cat. No.	Material	Description	Packaging
61080	PP	5 x 16 (80 Holes)	25

### 4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.

### 5. Product Photographs



SPL-TDS-MYURA


*For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.*  
 © SPL Life Sciences Co., Ltd., Version 2, July 28<sup>th</sup>, 2016

## 6. Specifications & Technical Drawing

### - Specifications

Dimensions		Material (Resin)	Packaging
Length:	228.00 mm	Polypropylene (PP)	▪ Pieces/Case: 25
Width:	63.40 mm		
Height:	27.55 mm		

### - Technical Drawing



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-MYURA

*For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.*  
© SPL Life Sciences Co., Ltd., Version 2, July 28<sup>th</sup>, 2016