

Incu Tissue

(#310070/310071/310072/310074)

1. Description

Incu Tissues are ideal for plant tissue cultures. As the sprout grows, another Incu Tissue can be connected on top as an optional frame (Cat. No. 310074).

2. Features

- ❖ Modular design
- ❖ Embossing on the closing flange of cover and the jar
- ❖ Three different materials (PS/PC/PP)
- ❖ Cat. No. 310070 can be assembled with Insect Breeding Box (Cat. No. 310075, 310076, 310077)
- ❖ Optical frame (Cat. No. 310074)

3. Product Range

Cat. No.	Material	Style	External Dimension w x l x h (mm)	Internal Dimension w x l x h (mm)	Sterile	Packaging
310070	PS	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	+	4 / 120
310071	PC	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	-	4 / 120
310072	PP	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	-	4 / 120
310074	PP	Frame	80.80 x 80.80 x 21.00	63.80 x 63.80 x 21.00	-	30 / 60

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-ICTI

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 21th, 2016


6. Specifications & Technical Drawing

(1) Square (#310070/310071/310072)

- Specifications

Dimensions	Material (Resin)			Packaging
	310070	310071	310072	
Width _{top} :	72.00 mm			<ul style="list-style-type: none"> • Pieces/Sleeve: 4 • Pieces/Case: 120
Length _{top} :	72.00 mm			
Width _{bottom} :	65.00 mm			
Length _{bottom} :	65.00 mm			
Height:	99.20 mm			
	Crystal-grade polystyrene (PS)	Polycarbonate (PC)	Polypropylene (PP)	

- Technical Drawing




(2) Frame (#310074)

- Specifications

Dimensions	Material (Resin)	Packaging
Width :	Polypropylene (PP)	<ul style="list-style-type: none"> • Pieces/Sleeve: 30 • Pieces/Case: 60
Length :		
Height:		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-ICTI

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© SPL Life Sciences Co., Ltd., Version 2, July 21st, 2016