

## Plant Culture Dish

(#310103/310100/310101/310200/10090/10091/10100)

### 1. Description

Plant Culture Dishes are used for sprout culture. Pre-sterilize, flat dishes provide excellent culture results.

### 2. Features

- ❖ Excellent air circulation design
- ❖ Excellent flatness
- ❖ Stable stacking
- ❖ Designed for reduced contamination
- ❖ Autoclavable (Cat. No. 310103, 310101)
- ❖ Gridded bottom (Cat. No. 310200)

### 3. Product Range

Cat. No.	Material	Dish Style d x h (mm)	Internal Dimension d x h (mm)	Sterile	Packaging
310103	PP	100.00 x 50.00	91.40 x 48.52	-	5 / 200
310100	PS	100.00 x 40.00	91.3 x 38.20	+	5 / 200
310101	PP	100.00 x 40.00	91.3 x 38.20	-	5 / 200
310200	PS	100.00 x 20.00	94.6 x 18.6	+	10 / 200
10090	PS	90.00 x 15.00	85.90 x 12.60	+	10 / 500
10091	PS	90.00 x 20.00	86.20 x 17.70	+	10 / 200
10100	PS	100.00 x 15.00	85.90 x 12.60	+	10 / 500

### 4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

### 5. Product Photographs



SPL-TDS-PCD

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.  
© SPL Life Sciences Co., Ltd., Version 2, July 21<sup>th</sup>, 2016


## 6. Specifications & Technical Drawing

### (1) Plant Culture Dish Ø100 (# 310103)

#### - Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 101.90 mm Height: 10.60 mm	Diameter: 102.00 mm Height: 50.00 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 5</li> <li>• Pieces/Case: 200</li> </ul>

#### - Technical Drawing




### (2) Plant Culture Dish Ø100 (# 310100/310101)

#### - Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 95.36 mm Height: 10.30 mm	Diameter: 91.35 mm Height: 40.00 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 5</li> <li>• Pieces/Case: 200</li> </ul>

#### - Technical Drawing




### (3) Plant Culture Dish Ø100 (# 310200)

#### - Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 96.57 mm Height: 9.00 mm	Diameter: 94.62 mm Height: 20.00 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>• Pieces/Sleeve: 10</li> <li>• Pieces/Case: 200</li> </ul>

#### - Technical Drawing




SPL-TDS-PCD

## (4) Plant Culture Dish Ø90 (# 10090)

## - Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 92.40 mm Height: 8.76 mm	Inner Diameter : 85.90 mm Inner Height : 12.60 mm Growth Area: 57.50 cm <sup>2</sup>	Crystal-grade polystyrene (PS)	<b>10090</b> Pieces/Sleeve : 10 Pieces/Case : 50

## - Technical Drawing




## (5) Plant Culture Dish Ø90 (# 10091)

## - Specifications

Dimensions		Material (Resin)	Packaging
Lid	Dish		
Diameter: 93.10 mm Height: 9.30 mm	Inner Diameter : 86.20 mm Inner Height : 17.70 mm Growth Area: 58.00 cm <sup>2</sup>	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>▪ Pieces/Sleeve: 10</li> <li>▪ Pieces/Case: 200</li> </ul>

## - Technical Drawing




SPL-TDS-PCD

## (6) Plant Culture Dish Ø100 (# 10100)

## - Specifications

Lid		Dish		Material (Resin)	Packaging
Dimensions					
Diameter:	102.20 mm	Inner Diameter :	96.40 mm	Crystal-grade polystyrene (PS)	<ul style="list-style-type: none"> <li>▪ Pieces/Sleeve: 10</li> <li>▪ Pieces/Case: 500</li> </ul>
Height:	8.85 mm	Inner Height :	13.75 mm		
		Growth Area:	72.30 cm <sup>2</sup>		

## - Technical Drawing



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-PCD

*For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.*  
© SPL Life Sciences Co., Ltd., Version 2, July 21<sup>th</sup>, 2016